
























Approved Denture Tooth Libraries

Approved Denture Tooth Libraries				Production Options						
	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds	 Monoblock for Try-In	 Printed base & Printed teeth	 Printed base & Pre-manufactured teeth	 Printed base & Milled teeth	 Milled base & Pre-manufactured teeth	 Milled base & Trimmed pre-manufactured teeth	 Milled base & Milled teeth
 Candolor										
Bonartic® II NFC+*	0	0	3	●	●	●	●	●	●	●
Condyliform® II NFC+*	0	0	3	●	●	●	●	●	●	●
PhysioStar® NFC+*	15	5	0	●	●	●	●	●	●	●
 Dentsply Sirona Inc.										
Genios® including 1to1 and 1to2 posterior setups	17	10	18	●		●		●	● (compatible with DS workflow)	
Portrait® IPN®	22	12	13	●		●		●	● (compatible with DS workflow)	
IPN 3D™ Portrait® Inspired (coming soon)										
 gebdi GmbH										
Tribos 501*	7	4	7	●	●	●	●	●	●	●
<p>* This tooth library contains open files that can be printed or milled. It can be selected whether they are produced as single teeth or as bridged teeth.</p>										

ADMIM-Full_denture_libraries-en, revision 2020-08-06

Approved Denture Tooth Libraries				Production Options						
	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds	 Monoblock for Try-In	 Printed base & Printed teeth	 Printed base & Pre-manufactured teeth	 Printed base & Milled teeth	 Milled base & Pre-manufactured teeth	 Milled base & Trimmed pre-manufactured teeth	 Milled base & Milled teeth
 KULZER <small>MITSUBI CHEMICALS GROUP</small> Kulzer GmbH										
Artic® digital	8	3	3	●		●		●	●	
Delara®	18	5	4	●		●		●	●	
Pala® Idealis (coming soon)										
Pala® Mondial (coming soon)										
 Myerson LLC.										
Kenson®*	18	9	9	●	●	●	●	●	●	●
 Nobilium®										
NOBILDENT™*	13	8	4	●	●	●	●	●	●	●
 Polident GmbH										
REF-LINE	7	1	6			●		●	●	
* This tooth library contains open files that can be printed or milled. It can be selected whether they are produced as single teeth or as bridged teeth.										

Approved Denture Tooth Libraries				Production Options						
	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds							
 Ruthinium® Group <i>We take care of your smile</i>										
Acry PLUS* including lingualized setups for all posterior moulds	28	12	10	●	●	●	●	●	●	●
Acry SMART	3	2	4			●		●	●	
 SHOFU® Inc.										
Veracia SA	9	4	4			●		●	●	
VITA VITA Zahnfabrik										
VIONIC VIGO® including buccal, crossbite, lingualized and Dr. Gerber setups	8	4	16	● (noised tooth shapes)		●		●	●	
VITAPAN EXCELL®	9	4	0	● (noised tooth shapes)		●		●	●	
VITAPAN LINGOFORM® including buccal, crossbite, lingualized and Dr. Gerber setups	0	0	16	● (noised tooth shapes)		●		●	●	
HUGE HUGE (coming soon)										
Sonning KF										
* This tooth library contains open files that can be printed or milled. It can be selected whether they are produced as single teeth or as bridged teeth.										

Approved Denture Tooth Libraries				Production Options						
	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds							
 Merz Dental GmbH (coming soon)										
artegral® Low resolution open data – Precise tooth pockets not possible										
 Unidesa Odi (coming soon)										
odident+										
* This tooth library contains open files that can be printed or milled. It can be selected whether they are produced as single teeth or as bridged teeth.										