
exocad GmbH

exoplan Supported Libraries,
solely with parts legally marketed in the U.S.
133 implant systems supported (2342 implants)

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
Anthogyr SAS (FR)						
Axiom® BL 2.8	3					
Axiom® BL PX	21					
Argon Dental Vertriebs GmbH & Co. KG (DE)						
K3Pro® Rapid	28	●				
K3Pro® Short	12	●		FG not for Ø 6.0 mm		
K3Pro® Sure	20	●				
Bicon®, LLC (US)						
MAX 2.5™ Integra-CP™	6	●		FG not for Ø 4.0 mm, length 8 mm		
Narrow® Integra-CP™	4	●		FG for Ø 3.0 mm		
SHORT® Integra-CP™	10	●		FG not for Ø 6.0 mm		
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available						

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
CAMLOG Biotechnologies GmbH (CH)						
CAMLOG® Screw-Line Promote® former-packaging	19		●			
CAMLOG® Screw-Line Promote® Plus Guided	11	●			● FG	
CAMLOG® Screw-Line Promote® Plus former-packaging	19		●			
CAMLOG® Screw-Line Promote® Plus snap-in	19		●			
CAMLOG® Screw-Line Promote® snap-in	19		●			
CONELOG® Screw-Line Promote® Plus Guided	14	●			● FG	
CONELOG® Screw-Line Promote® Plus former-packaging	19		●			
CONELOG® Screw-Line Promote® Plus snap-in	19		●			
iSy®	12					
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available						

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
CSM Implant Co., Ltd. (KR)						
Submerged III	21					
Submerged	20					
Dentsply Sirona Inc. (US)						
Ankylos® C/X	22	●		FG not for Ø 3.5 mm + Ø 4.5 mm, length 6.6 mm + 17 mm; FG not for Ø 5.5 mm + Ø 7.0 mm	● FG	
Astra Tech OsseoSpeed™ EV	44					
Astra Tech OsseoSpeed™ EV Profile	24					
Frialit®	17					
XiVE® S	25	●		FG not for Ø 3.4 mm, Ø 3.8 mm + Ø 4.5 mm, length 18 mm; FG not for Ø 5.5 mm	● FG	
XiVE® TG	14					
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available						

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
Leone® S.p.A. (IT)						
2.9	3					
Implant	12					
Max Stability	8					
Short 6.5	1					
medical instinct® (DE)						
BoneTrust® plus cone	16					
BoneTrust® plus hex	16					
MIS Implants Technologies Ltd. (IL)						
C1	19	●	● UNI	FG for implant lengths 8 mm - 13 mm		
M4	30					
SEVEN	30	●	● UNI	FG for implant lengths 6 mm - 13 mm; FG not for Ø 6.0 mm		
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available						

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
V3	19	●	● UNI	FG for implant lengths 8 mm - 13 mm		
OSSTEM IMPLANT CO., LTD (KR)						
SS Implant System	26					
TS Implant System	35	●		FG for implant lengths 7 mm - 13 mm (FG not for 7 mm extra short)		
US Implant System	18					
SIC invent AG (CH)						
SICace®	20					
SICmax® onepiece	2					
SICmax®	20					
SICvantage® max	24					
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available						

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
S.I.N. - Sistema de Implante Nacional S.A. (BR)						
REVOLUTION®	4					
STRONG SW® CM	20					
STRONG SW® HE	17					
STRONG SW® HI	15					
TRYON®	23					
UNITITE COMPACT®	9					
UNITITE PRIME®	15					
UNITITE SLIM®	3					
Southern Implants (ZA)						
Deep Conical Cylindrical	18					
Deep Conical Tapered	18					
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available						

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
External Hex - Cylindrical	34		● FG	PT for implant lengths 8.5 mm - 15 mm; PT not for Ø 4.0 mm		
External Hex - Cylindrical-internal drive	21					
External Hex - Platform Switch	4					
External Hex - Tapered	32	●	● FG	FG for Ø 3.25 mm, Ø 4.0 mm + Ø 4.7 mm, lengths 8.5 mm - 15 mm; PT for implant lengths 8.5 mm - 15 mm		
External Hex - Tapered-Internal Drive	14					
IT-Cylindrical	19					
IT-Tapered	14					
Provata - Tapered	12	●		FG not for length 18 mm		
Tri-Nex - Cylindrical	15					
Tri-Nex - Tapered	20					

FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
Institut Straumann® AG (CH)						
Bone Level Roxolid® SLA® Loxim®	12					
Bone Level Roxolid® SLActive® Guided	9	●	● FG		● FG	● FG
Bone Level Roxolid® SLActive® Loxim®	12					
Bone Level Tapered Roxolid® SLA® Loxim®	15					
Bone Level Tapered Roxolid® SLActive® Guided	9	●	● FG		● FG	● FG
Bone Level Tapered Roxolid® SLActive® Loxim®	18					
Bone Level Tapered Titanium SLA® Loxim®	15					
Bone Level Titanium SLA® Loxim®	12					
Bone Level Titanium SLA®	12					
PURE Ceramic Implant	4					
PURE Ceramic Monotype Implant	16					

FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
Standard Plus Roxolid® SLA® Loxim®	22					
Standard Plus Roxolid® SLActive® Guided	9	●	● FG		● FG	● FG
Standard Plus Roxolid® SLActive® Loxim®	22					
Standard Plus Titanium SLA® Loxim®	18					
Standard Plus Titanium SLA®	18					
Standard Roxolid® SLA® Loxim®	20					
Standard Roxolid® SLActive® Guided	6	●	● FG		● FG	● FG
Standard Roxolid® SLActive® Loxim®	20					
Standard Titanium SLA® Loxim®	20					
Standard Titanium SLA®	20					
Tapered Effect Roxolid® SLA® Loxim®	11					
Tapered Effect Roxolid® SLActive® Guided	3	●	● FG		● FG	● FG

FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available

ADMIM-exoplan-libraries_US-en, revision 2019-08-23

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
Tapered Effect Roxolid® SLActive® Loxim®	11					
Tapered Effect Titanium SLA® Loxim®	11					
Tapered Effect Titanium SLA®	11					
TAG Dental (IL)						
Axis	22	●		FG for Ø 3.3 mm, Ø 3.75 mm + Ø 4.2 mm, lengths 8 mm - 13 mm		
Massif	18					
Thommen Medical AG (CH)						
CONTACT implant MC INICELL®	19					
CONTACT implant®	21					
CONTACT implant RC INICELL®	21					
ELEMENT implant LC INICELL®	11					
ELEMENT implant MC INICELL®	24					
FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available						

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
ELEMENT implant®	30					
ELEMENT implant RC INICELL®	34					
TRI® Dental Implant System Int. AG (CH)						
TRI®-Narrow	3					
TRI®-Octa	14					
TRI®-Vent	17					
Uris Inc. (US)						
OMNI System	23	●				
Zimmer Biomet (US)						
Certain® 3i T3® for Parallel Walled Navigator®	22					
Certain® 3i T3® with DCD® for Parallel Walled Navigator®	22					
Certain® Full OSSEOTITE® for Parallel Walled Navigator®	20					

FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
Certain® NanoTite® for Parallel Walled Navigator®	22					
Certain® OSSEOTITE® for Parallel Walled Navigator®	23					
Certain® OSSEOTITE®2 for Parallel Walled Navigator®	23					
Certain® PREVAIL® 3i T3® for Parallel Walled Navigator®	15					
Certain® PREVAIL® 3i T3® with DCD® for Parallel Walled Navigator®	15					
Certain® PREVAIL® NanoTite® for Parallel Walled Navigator®	10					
Certain® PREVAIL® OSSEOTITE® for Parallel Walled Navigator®	10					
Certain® PREVAIL® OSSEOTITE®2 for Parallel Walled Navigator®	15					
Certain® Tapered 3i T3® for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm		
Certain® Tapered 3i T3® with DCD® for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm		
Certain® Tapered CP4 TAPERED for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm		
Certain® Tapered Full OSSEOTITE® for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm		

FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available

ADMIM-exoplan-libraries_US-en, revision 2019-08-23

Available Implant Systems		Available Original Sleeves			Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco
Certain® Tapered NanoTite® for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm		
Certain® Tapered OSSEOTITE® for Tapered Navigator®	20	●		FG for Ø 3.25 mm, Ø 4.0 mm + Ø 5.0 mm		
Certain® Tapered PREVAIL® 3i T3® for Tapered Navigator®	15	●		FG for Ø 4.0 mm + Ø 5.0 mm		
Certain® Tapered PREVAIL® 3i T3® with DCD® for Tapered Navigator®	15	●		FG for Ø 4.0 mm + Ø 5.0 mm		
Certain® Tapered PREVAIL® CP4 TAPERED for Tapered Navigator®	15	●		FG for Ø 4.0 mm + Ø 5.0 mm		
Certain® Tapered PREVAIL® NanoTite® for Tapered Navigator®	15	●		FG for Ø 4.0 mm + Ø 5.0 mm		
External Hex 3i T3®	34					
External Hex 3i T3® with DCD®	30					
External Hex Full OSSEOTITE®	29					
External Hex NanoTite®	27					
External Hex OSSEOTITE®	34					
External Hex OSSEOTITE®2	30					

FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available

Available Implant Systems		Available Original Sleeves				Compatible Sleeves for FG	
Implant system	Included implants	Fully guided	Pilot	Comments	nt-trading	Steco	
External Hex Tapered 3i T3®	20						
External Hex Tapered 3i T3® with DCD®	20						
External Hex Tapered Full OSSEOTITE®	20						
External Hex Tapered NanoTite®	20						
External Hex Tapered OSSEOTITE®	20						

FG = Fully guided workflow available, UNI = Universal compatibility, PT = Compatible pilot sleeves available

Sleeve Manufacturer	Available Sleeves		
	Fully guided	Universal	Pilot
3D Diagnostix, Inc. (US)		Compatible with: all implants	Compatible with: all implants
Argon Dental Vertriebs GmbH & Co. KG (DE)	Compatible with: Argon Dental Vertriebs GmbH & Co. KG		
Bicon®, LLC (US)	Compatible with: Bicon®, LLC		
CAMLOG Biotechnologies AG (CH)	Compatible with: CAMLOG Biotechnologies GmbH		Compatible with: CAMLOG Biotechnologies GmbH
C-Tech Implant S.R.L. (IT)	Compatible with: C-Tech Implant S.R.L.		
Dentsply Sirona Inc. (US)	Compatible with: Dentsply Sirona Inc.		
Medentika® GmbH (DE)	Compatible with: Medentika® GmbH		
MIS Implants Technologies Ltd. (IL)	Compatible with: MIS Implants Technologies Ltd.		Compatible with: MIS Implants Technologies Ltd.
nt-trading & GmbH Co. KG (DE)	Compatible with: CAMLOG Biotechnologies GmbH, Dentsply Sirona Inc., Institut Straumann® AG		Compatible with: all implants, Institut Straumann® AG
OSSTEM IMPLANT CO., LTD (KR)	Compatible with: OSSTEM IMPLANT CO., LTD		
Southern Implants (ZA)	Compatible with: Southern Implants		Compatible with: Southern Implants
steco-system-technik GmbH (DE)	Compatible with: Institut Straumann® AG	Compatible with: all implants	Compatible with: all implants
Institut Straumann® AG (CH)	Compatible with: Institut Straumann® AG		Compatible with: Institut Straumann® AG
TAG Dental (IL)	Compatible with: TAG Dental		
Uris Inc. (US)	Compatible with: Uris Inc.		
Zimmer Biomet (US)	Compatible with: Zimmer Biomet		